



AGH UNIVERSITY OF SCIENCE
AND TECHNOLOGY

DETERMINATION OF RADIUM ISOTOPES OF VERY LOW CONCENTRATIONS IN WATER SAMPLES BY LSC TECHNIQUES

Nguyen Dinh Chau, Gargul Magdalena

Faculty of Physics and Applied Computer Science
Department of Applied Nuclear Physics

Niš, Serbia - 29.05.2014



Outlines of presentation

- Radiological criteria for drinking water
- Guidance levels for natural radionuclides in drinking water
- Critical Level and Low Limit Detection
- Chemical and measurement procedures for radium isotopes determination
- Factors controlling the background level in the LSC method
- Influence of BaCl_2 on the measured background samples
- Influence of BaCl_2 on chemical efficiency
- Test of the method
- Conclusions

Radiological criteria for drinking water

- According to WHO Guidelines the annual committed effective dose resulting from absorb of all radionuclides in drinking water should be below 0.1 mSv/year
- According to Guideline, the max level (GL) for radionuclides in drinking water can be calculated as follow:

$$GL = \frac{IDC}{h_{ing}(ag) \cdot V(ag)}$$

where:

GL – maximum level of radionuclide in drinking water [mBq/L],

IDC – annual permissible committed dose; $IDC = 0.1 mSv$,

$h_{ing}(ag)$ – dose coefficient for ingestion by person in the ag – age group [mSv/mBq],

$V(ag)$ – annual ingested volume of drinking waters

$$\sum_i \frac{C_i}{GL_i} \leq 1$$

Sum of all ratios of the real concentration of i-isotope to its GL should be lower than 1.

Guidance levels for some selected natural radionuclides in drinking water [*mBq/L*]



Isotope	Age group					
	< 1 y	1-2 y	2-7 y	7-12 y	12-17 y	> 17 y
U-238	403	1142	1712	2014	2045	3044
U-234	370	1054	1557	1851	1851	2795
Th-232	30	304	391	472	548	595
Ra-226	29	143	221	171	91	489
Ra-228	5	24	40	35	25,6	198
Ra-224	51	208	391	527	685	2107
Pb-210	1.6	38	62	72	72	198
Po-210	5	16	31	53	86	114

Critical level and limit of detection

Follow Lloyd Currie:

Critical level: $L_C = k\alpha \cdot \sigma_B$

Detection limit: $L_D = L_C + k\beta \cdot \sigma_D$



Follow ISO 11929:2010:

$$g(L_i) = L_i$$


$$h(L_i) = L_i + k \cdot \sigma_i$$



The LLDs requirement of Council Directive for analyzing methods and LLDs of routine methods

Isotope	COUNCIL DIRECTIVE 2013/51/EURATOM	Routine methods
U-238	0.02 Bq/L	0.005 Bq/L
U-234	0.02 Bq/L	0.005 Bq/L
Ra-226	0.04 Bq/L	0.01 Bq/L
Ra-228	0.02 Bq/L	0.03 – 0.8 Bq/L
Pb-210	0.02 Bq/L	0.03 – 0.8 Bq/L
Po-210	0.01 Bq/L	0.01 Bq/L

Chemical procedure for radium in water sample



The chemical agents used in chemical procedure:

- HNO_3
- EDTA
- $\text{Pb}(\text{NO}_3)_2$
- BaCl_2
- H_2SO_4
- CH_3COOH

In rououtine analyzing, the background sample is to be prepared together with analyzed samples using distilled water.



LSC spectrometers

Guardian 1414 α/β Liquid Scintillation counter



Quantulus



RAD2014 - Niš, Serbia




Factors controlling the background level in LSC

- **The electronic noise;**
- **The cosmic and radiation ray from surrounding;**
- **Chemical Reagents;**
- **Scintillation cocktail and glass vial.**

The influences of the chemical reagent BaCl_2 were carried out; the blank samples were prepared using 0.5 L of distilled water with different amounts of BaCl_2 , but the amounts of the other chemical reagents were the same as described in procedure.

Measured count rates of the α and β canals for blank samples with different masses of BaCl_2 -0.1M



Concentration of Ra-226 in chemical reagent BaCl_2 – 0.1 M

Volume of BaCl_2 -0.1 M	18 mL	16 mL	14 mL	12 mL
Average α nett count rates (cpm)	7.2	6.6	6.2	5.3
Stand. Uncert.	1.2	0.9	1.3	0.8

Based on the data in the above table, the concentration of Ra-226, amounts to **84 ± 8 Bq/kg** of BaCl_2 . This value was confirmed through gamma spectrometry.





Critical level (L_c)

$$L_c = 1.65 \cdot \sigma_0 *$$

BaCl ₂ -0.1 M	18 mL	16 mL	14 mL	12 mL	10 mL
I α (cpm)	7.2 \pm 1.2	6.6 \pm 0.9	6.2 \pm 1.3	5.3 \pm 0.8	5.0 \pm 0.8
I β (cpm)	15.4 \pm 3.1	15.8 \pm 2.9	14.2 \pm 2.8	12.2 \pm 2.2	12.0 \pm 2.5
Ra-226 (mBq)	8	6	8	6	5
Ra-228 (mBq)	33	30	28	20	25


* by Lloyd Currie

Chemical efficiency as a function of the added amounts BaCl_2 0.1 M



Optimal amount of BaCl_2 – 0.1 M

$$F = \frac{\text{signal}^2}{\text{background}}$$



Test of the method

Sample	Spiked contents of radium isotopes [mBq]		Determined [mBq]	
	$^{226}\text{Ra} \pm \sigma$	$^{228}\text{Ra} \pm \sigma$	$^{226}\text{Ra} \pm \sigma$	$^{228}\text{Ra} \pm \sigma$
S1	0	23 ± 2	nd* (≤ 5.0)	17 ± 3
S2	51 ± 2	0	48 ± 2	nd (≤ 20)
S3	53 ± 2	52 ± 3	48 ± 2	50 ± 8



Conclusions

- The current reagent of BaCl_2 is Ra-226 contaminated with 84 Bq/kg;
- Critical level for determination of Ra-226 and Ra-228 depend on the used amount of BaCl_2 ;
- 12 mL BaCl_2 0.1 M solution is the optimal used amount to preparing;
- The LSC method is still to be improved, unless the method only can be used for Ra-226 analyzing;
- For high precise determination of Ra-228 we should combine the gamma and LSC methods.



RAD2014 - Niš, Serbia